



DENTAL FLOSS AID INSTRUCTION MANUAL



Flossing should not feel like a workout for your hands.

This dental flosser aid is designed to make flossing easier and more manageable, especially if grip, dexterity, or hand strength are a challenge. The handle gives you something solid to hold onto, while the open fork design keeps floss taut and easy to guide between teeth.

It's a simple tool that can make daily oral care feel more achievable, whether that's because of arthritis, tremors, chronic pain, muscle weakness, or just wanting a quicker, easier way to floss. It's also a great option for kids who are learning to floss independently.

Lightweight and portable, it's easy to keep at home or take with you when you're out and about.

Great for:

- Making daily flossing easier and more accessible
- Reduced hand strength or grip fatigue
- Dexterity challenges or tremors
- Arthritis, joint pain, or stiffness in the hands
- Chronic pain or hand fatigue
- Neurological conditions that affect fine motor control
- Limited finger movement or coordination
- Oral care after dental work or sensitivity
- Children learning independent flossing skills
- Older adults who find flossing difficult
- Portable, on-the-go oral care
- Simplifying daily routines at home
- Anyone who finds traditional flossing awkward, tiring, or frustrating



DENTAL FLOSS AID INSTRUCTION MANUAL

Features:

- Holds dental floss securely and keeps it taut
- Easy-to-hold handle for better control
- Wide fork design to reach between teeth
- Reusable and washable
- Lightweight and compact
- Available in green or blue

Image Description: Plastic dental flosser aids shown in blue and green. The flossers have a flat handle with a circular centre dial and a forked top designed to hold floss securely across the open ends. Images show the flosser from multiple angles and in hand to demonstrate size and grip.

Details:

Material: PP Plastic

Size: 12cm (L) x 3cm (W)

Weight: 30 grams

Package Contains: 1x Dental Flosser Aid

Age Guidance: Suitable for ages 3 and up

How to Use:

1. Cut a length of dental floss
2. Wrap one end around the centre dial a few times
3. Bring the floss up and over the arms of the flosser
4. Wrap the other end around the dial to keep the floss tight
5. Floss between teeth as usual
6. Remove or slide the floss out when finished

Tip: Use gentle pressure and avoid forcing the flosser

Care Instructions:

- Rinse with warm water after use
- Wash with mild soap as needed
- Air dry before storing
- Do not use boiling water