



DENTAL WATER FLOSSER + IRRIGATOR, 6 PIECE

INSTRUCTION MANUAL



Flossing does not have to be fiddly or uncomfortable. This electric dental flosser uses a gentle stream of water to do the hard work for you, helping clean between teeth and along the gum line without traditional string floss. Simply fill the water tank, choose your attachment, and select the mode that feels right for you. The water based flow helps remove food debris and plaque while being gentler on gums and less abrasive than regular flossing. It is especially helpful if you struggle with dexterity, jaw movement, or reaching awkward angles in your mouth. With three pressure modes and a built in safety shut off, this flosser supports daily dental care in a way that feels more manageable, comfortable, and independent.

Great for:

- Braces, retainers, or dentures where food build up is common
- People with bleeding disorders or low platelet levels
- Sensory sensitivities around traditional flossing
- Limited grip strength, dexterity, or range of motion
- Wrist, arm, or hand pain
- TMJ or reduced jaw opening
- Neck or spinal limitations, including fusions
- Disabled people who need easier access to dental hygiene
- Anyone wanting a gentler, more thorough flossing routine

Features:

- Water based flossing for gentle, effective cleaning
- Three modes: soft, strong, and pulse
- Auto shut off after 2 minutes for safety
- Fully waterproof design
- Attachments store neatly inside the water cup
- Lightweight and travel friendly

Image Description: An electric oral flosser shown with multiple interchangeable nozzle attachments. The main unit is white with simple control buttons. Close up images show the water tank, angled tips, and charging cable.



DENTAL WATER FLOSSER + IRRIGATOR, 6 PIECE

INSTRUCTION MANUAL

r

Details:

Material: Food Grade ABS Plastic

Size: 28cm (L) x 6cm (Diameter)

Water Capacity: 220ml

Weight: 250 grams

Package includes: 1x Electric Dental Oral Flosser, 4x dental flossing attachments, + Instructions Manual

Age Guidance: Suitable for teens and adults. Not recommended for young children without supervision.

How to Use:

1. Fill the water tank and secure the lid
2. Choose and attach your preferred nozzle
3. Select your desired mode
4. Gently guide the water stream between teeth
5. Empty and rinse after use

Care Instructions:

- Rinse attachments after each use
- Empty the water tank when finished
- Wipe the exterior with a damp cloth
- Allow to dry fully before storing