

ADJUSTABLE CATHETER BAG HOLDER



Description:

Having a catheter bag shouldn't mean feeling exposed, restricted, or unsure about where to put it. This adjustable catheter bag holder is designed to support comfort, discretion, and confidence in everyday life. It keeps your catheter bag upright and secure while giving you flexible options for how and where you carry it, whether that's over your shoulder, attached to a wheelchair, or secured to a bed frame.

With a dual-sided design, one side offers a water-resistant layer for privacy and protection, while the other features breathable mesh to allow airflow and help reduce moisture build-up. Lightweight and easy to adjust, it's a practical solution that supports independence at home, in hospital, or when you're out and about.

Great for:

-  Supporting everyday catheter use with confidence
-  Maintaining privacy and discretion
-  Attaching to wheelchairs, beds, or mobility equipment
-  Hospital stays or recovery at home
-  Travel and day-to-day outings
-  Supporting independence and hygiene

Features:

-  Dual-sided design with water-resistant fabric and breathable mesh
-  Adjustable straps with strong clips for secure attachment
-  Mesh panel helps with airflow and moisture control
-  Lightweight and easy to carry or reposition
-  Multiple wearing and hanging options
-  Keeps catheter bags upright to support hygiene and reduce leaks

Details:

Material: Canvas

Size: *Please see below...*

Bag: 27cm (L) x 19cm (W)

Strap: Adjustable up to 110cm (L)

Weight: 83 grams

Package Contains: 1x Catheter Bag Holder

Age Guidance:

Please Note: colour may vary slightly...

How to Use:

-  Place your catheter bag securely inside the holder
-  Choose the side you prefer facing out, water-resistant or mesh
-  Adjust the strap length to suit your comfort
-  Clip onto a wheelchair, bed frame, or wear over the shoulder
-  Ensure the bag hangs upright and comfortably during use
-  Regularly check the mesh, straps, stitching and buckles for wear or damage. Replace if damaged or no longer holding the drainage bag securely.
-  Keep the drainage bag positioned below bladder level and make sure the holder or straps do not restrict circulation, kink the tubing or pull on the catheter.

Care Instructions:

-  Wipe clean with a damp cloth as needed
-  Hand wash gently if required
-  Allow to air dry fully before reuse