



UMBRELLA CLAMP HOLDER INSTRUCTION MANUAL



Rain or shine, this umbrella holder has you covered, literally. It clips onto wheelchairs, walkers, mobility scooters, strollers, or bikes, so you're not trying to juggle an umbrella while you're moving. The clamp holds steady, and the adjustable arm lets you shift your umbrella into the position you want for shade or shelter.

Set it once, and get on with what you're doing.

Great For:

- Protecting yourself from harsh sun on hot days
- Staying dry when the rain starts unexpectedly
- Wheelchair users who need both hands free
- Walkers who want shade or shelter on the go
- Mobility scooter riders
- Parents pushing strollers
- Biking with light rain or sun cover
- Outdoor markets, events, and festivals
- Picnics, sports days, or community outings
- Shopping trips when you don't want to fuss with holding an umbrella
- Travel and sightseeing where weather can change quickly

Features:

- Strong clamp grips securely onto round frames or bars
- Adjustable arm lets you change angle and position
- Fits a range of mobility devices, bikes, and strollers
- Lightweight design that doesn't add bulk
- Simple to attach and remove when needed



UMBRELLA CLAMP HOLDER INSTRUCTION MANUAL

Image Description: A hands-free umbrella holder with a sturdy black clamp that attaches securely to wheelchairs, scooters, strollers, or bikes. The flexible, adjustable arm extends upward to hold an umbrella, keeping it stable while allowing the user to stay shaded or dry without needing to hold it.

Details:

Material: Stainless Steel + ABS Plastic

Size: 36cm (L) x 3cm (diameter tube) approx

Weight: less than 250 grams

Package includes: 1x Folding Umbrella Clamp

Age Guidance: Suitable for teens and adults. Not recommended for young children due to clamp mechanism and moving parts. Adult supervision recommended where needed. Not a toy.

How to Use:

1. Open the clamp and position it around a round frame or bar
2. Tighten the clamp until secure
3. Insert your umbrella into the holder
4. Adjust the arm to your preferred angle

Care Instructions:

- Wipe clean with a damp cloth
- Dry fully after exposure to rain
- Do not overload or force beyond natural bend
- Store in a dry place when not in use

Safety Information:

- Ensure clamp is securely fastened before use
- Check stability before moving
- Not suitable for strong winds or heavy-duty conditions
- Use caution when adjusting near face or eyes